



Esco Sample Protection and Remote Monitoring Solutions *Sample monitoring made easy*

- **Voyager®** (*Remote Monitoring, Programming, Datalogging Software*)
- **iRMAP™** (*Intelligent Remote Monitoring Application Protocol*)
- **Voyager® Box** (*Web Monitoring and Datalogging System*)
- **PROtect** (*A Completely Independent and Redundant Sample Monitoring System*)
- **SampleTracker** (*Biosample Inventory Management System*)

Esco Sample Protection and Remote Monitoring Solutions

Table of Contents

Corporate Profile	3
Products and Applications	4
Introduction	5
Voyager® Remote Monitoring, Programming, Datalogging Software	6
Voyager® Key Features	7
Voyager® Interface with Esco Equipment	8
Voyager® Installation and System Requirements	9
iRMAP™ Intelligent Remote Monitoring Applications Protocol	10
iRMAP™ Key Features	11
Voyager® Box	13
Voyager® Box Features	14
Voyager® Box Installation and System Requirements	15
PROtect A Complete Independent and Redundant Sample Monitoring System	16
PROtect Key Features	17
PROtect System Components and Requirements	19
PROtect General Specification	20
PROtect Options and Accessories	21
SampleTracker	24
SampleTracker Key Features	25
Esco After sales Service	27
Esco Global Network	28





Welcome to Esco

Esco's Vision is to provide enabling technologies for scientific discoveries to make human lives healthier and safer.



Esco represents innovation and forward-thinking designs, which are all coupled with the highest standard quality since 1978. The Esco Group of Companies remains dedicated in delivering innovative solutions for the clinical, life sciences, research, industrial, laboratory, pharmaceutical and IVF community. With the most extensive product line in the industry, our products have passed a number of international standards and certifications. Esco operates under ISO 9001, ISO 14001 and ISO 13485.

Availability and Accessibility. Headquartered in Singapore, manufacturing facilities are located in Asia and Europe. R&D is conducted worldwide spanning the US, Europe and Asia. Sales, services and marketing subsidiaries are located in 12 major markets including the US, UK, Japan, China and India. Our regional distribution centers are located in China, UK, India, Malaysia, Philippines, Singapore, South Africa, South Korea and United States of America. Because of our worldwide presence, you can be sure that Esco is within your reach.

High Quality, Reliable and Dependable. Our customers are confident that only with the best quality, reliable, and dependable products, can they be sure of the accuracy of their research and procedures. Cross functional teams from Esco Production, R&D, Quality Assurance and Senior Management, are regularly assembled to review and implement areas for improvement.

Esco Cares for Your Safety. Esco focuses on providing safety not just for your samples but also for users.

Esco Cares for Your Comfort. Comfort of our users is ensured by building ergonomic designs and by reducing the noise levels of the units.

Esco Cares for the Environment. One in every four of Esco's employees is involved in Research and Development and a number of these evaluate new components and/or designs to produce energy efficient equipment. Being GREEN is more than just modifying the parts we use to produce a new energy efficient technology, it also embodies the every aspect of the company.

Customer Service and Support. Our service does not stop once purchase has been made. Esco gives on-time customer service and offers end-user seminars, service training, preventive maintenance, provides educational materials and informative videos.

As Esco takes the opportunity to respond to the world's needs, we aim not just to contribute in the advancement of scientific discoveries but also in making the world a safer, healthier and a better place to live in.

Products and Application

Life Sciences Laboratory Equipment

Sample Preparation

- Class I Biological Safety Cabinets
- Class II Type A2 Biological Safety Cabinets
- Class II Type B1 Biological Safety Cabinets
- Class II Type B2 Biological Safety Cabinets
- Class III Biological Safety Cabinets
- Horizontal Laminar Flow Clean Benches
- Vertical Laminar Flow Clean Benches
- Laboratory Animal Research Workstations
- Laboratory Centrifuges

Sample Cultivation

- CO₂ Incubators, Direct Heat Air-Jacketed
- CO₂ Incubators with Cooling System
- CO₂ Incubators with Stainless Steel Exterior
- CO₂ Incubators (Water-Jacketed)
- Laboratory Shakers

Sample Analysis

PCR Thermal Cyclers

- Conventional Thermal Cyclers
- Real-time PCR Systems

PCR Sample Handling

- Microplate Shakers
- PCR Cabinets

Sample Storage & Sample Protection Solutions

- Ultra-low Temperature Freezers
- Lab Refrigerators and Freezers
- Sample Database Management Software
- Intelligent Remote Monitoring Application Protocol
- Remote Monitoring, Datalogging, Programming Software
- Wireless Monitoring System

Chemical Research

- Ductless Fume Hoods
- Laboratory Fume Hoods
- Fume Hood Airflow Monitors
- Exhaust Blowers
- Powder Weighing Balance Enclosures

General Equipment

Laboratory Thermostatic Products

- Laboratory Oven
- Laboratory Incubator
- Refrigerated Incubator
- Natural Convection Incubator

Forensic Sciences

- Evidence Drying Cabinet

Medical / IVF Equipment

Controlled Embryo Handling

- IVF Workstation
 - Stereo Zoom Microscope
- Anti-vibration Table

Safe Embryo Culture

- Benchtop Multi-room Embryo Incubators
- CO₂ Incubators

Innovative Time Lapse Imaging

- Time-Lapse Embryo Incubator

Accurate Quality Control

- CO₂ / O₂ Temperature Validation Unit

Healthcare

Esco Pharma Products

Airflow Containment Products

- Pharmacon® Downflow Booths
- Ceiling Laminar Airflow Units
- Laminar Flow Horizontal Trolley
- Laminar Flow Vertical Trolley
- Enterprise® Laminar Flow Straddle Units
- Garment Storage Cabinet
- Cytotoxic Safety Cabinets

Isolation Containment

- Aseptic Containment Isolator (ACTI)
- Weighing and Dispensing Containment Isolator (WDCI)
- General Processing Platform Isolator (GPPI)
- Containment Barrier Isolator (CBI)
- Turbulent Flow Aseptic (Grade A) Isolator (TFAI)
- Isoclean® Healthcare Platform Isolator (HPI)
- Streamline® Compounding Isolators (SCI)
- Technetium Dispensing Isolators
- Blood Cell Labeling Isolators
- Open and Closed Restricted Barrier Access Systems (RABS)

Cross Contamination Facility Integrated Barrier

- Cleanroom Air Showers
- Infinity® Air Shower Pass Box
- Infinity® Cleanroom Transfer Hatch
- Infinity® Pass Boxes
- Soft Capsule® Soft Wall Cleanroom
- Dynamic Passboxes and Dynamic Floor Label Hatches
- Esco BioPass™ Pass Through

Ventilation Containment

- Ventilated Balance Enclosure
- Extraction Hoods/Enclosures
- Local Exhaust Ventilation Systems

VacciXcell Products

Bioreactors and Fermenters

- CelCradle™
- TideCell®
- StirCradle™
- StirCradle™-Pro
- VacciXcell™ Hybrid bioreactor

Cell Culture Monitoring, Media and Consumables

- Super Plus™
- Plus™ Vero
- Plus™ MDCK
- Plus™ MDCK II
- BioNOC™ II macrocarriers
- GlucCell® Glucose Monitoring System
- CVD Kit

Filling Line Equipment

- AseptiCell®
- Traditional Filling Line

Integrated Solutions

- Cell Processing Isolator
- Cell Processing Center

TaPestle Rx Products and Services

PRODUCTS

Pharmacy Automation and Compounding Supply

- Compounding Pharmacy Isolators (SCI, HPI, CBI, GPPI)
- Safety Cabinets and Enclosures (Class II BSC, VBE, LFC)
- Radiopharmacy Hoods and Isolators
- Automated IV Compounding System*
- Aseptic Filling Systems

**Southeast Asian Markets only*

Healthcare and Laboratory Construction Components

- Prefabricated Walls (Airecell®)
- Prefabricated Containerized Facility (Prefab™)
- Series Ceiling Systems
- Hygienic/Hermetic Door Systems
- Surgical Scrub Sinks
- Vinyl Tiles and Epoxy
- Laboratory Fit-outs
 - Worktops
 - Frames
 - Specialty Storage cabinets
 - Service Spines & Reagent Shelving

SERVICES

- Conceptualization
- Planning
- Procurement
- Installation


FACILITY DESIGNS

- Process Architecture
- Biocontainment/Biosafety
- Pharmacy Compounding/Nuclear Medicine
- Cleanroom, Vaccine and Cell Processing
- Laboratory
- Containerized Facility
- ART/IVF
- Cold Chain

Every precaution has to be taken to protect a sample's integrity.

Be assured with the guaranteed quality of performance and service Esco can provide.

Esco's Sample Protection and Remote Monitoring Solutions provides an extensive solution to prevent loss and spoilage of unique samples and products.



Voyager®

Remote Monitoring, Programming, Datalogging Software

Esco Voyager® is a PC-based software package developed for the remote monitoring, datalogging, and programming/ device configuration of Esco controlled environment laboratory equipment. The equipment has the capability to push updates to Esco equipment connected to software. In any event of equipment failure, the software is able to send out remote alarm notification through e-mail or SMS with a modem configured. It is also compliant to 21 CFR Part 11 standard for traceability.

- Ability to generate reports and graph of different device parameters

- Protect the samples by providing alarms when exceeding a user-defined parameter limit via automatic email alerts

- Ability to save and report log data in various formats

- Easy PC interface using USB, LAN cables or via Wi-Fi

- Allows all Esco thermostatic products to be connected and monitored for proper operation from a central location



Compatible Esco Equipment (Clockwise from top):

- Lexicon® Ultra-low Temperature Freezer (UUS)
- Isotherm® Refrigerated Incubator (IFC)
- Isotherm® Forced Convection Incubator (IFA)
- CelCulture® CO₂ Incubator (CCL)

KEY FEATURES




Remote Monitoring

- Automatic, continuous monitoring of device parameters (Temperature, CO₂, RH, etc.) to validate instrument performance and assay criteria
- Automatic e-mailing of log data at user-defined intervals



Remote Programming and Equipment Configuration

- Development of programs using a graphical interface, for subsequent download to device memory
- Voyager® software can download user programs to multiple Esco instruments easily insuring consistent setup on multiple devices.
- Configure devices remotely (measurement interval, temperature set point, low temperature alarm point and high temperature alarm point)
- This window also allows the user to choose the method of when and which data can be viewed (spreadsheet, email, web page).



Datalogging and Graph Functions

- Easy viewing and graphing multiple device parameters remotely in real-time
- Saving and exporting log data in various formats from multiple instruments
- Download of data logs from equipment memory (only for devices with built-in data logging memory, i.e. CO₂ Incubators)



Alarm

- Alarm is activated if a user-defined parameter limit is exceeded.
- Automatic email alerts
- Documentation of all alarms with time and date stamp

Voyager® interfaces with individual Esco equipment over various available communication ports with the aid of EscoBUS communications protocol.

Compatible Equipment

- Lexicon II® – Ultra-Low Temperature Freezer (Silver, Gold, and Platinum)
- CelCulture® – CO₂ Incubator (CCL)
- CelMate® – CO₂ Incubator (CLM)
- Isotherm® – Forced Convection Laboratory Oven
- Isotherm® – Natural Convection Laboratory Incubator (INA)
- Isotherm® – Forced Convection Laboratory Incubator (IFA)
- Isotherm® – Refrigerated Incubator (IFC)



Installation Requirements



Compatible Equipment



Ethernet Cable



RS 485 Converter



Voyager® Software
installed in PC or Laptop

- Industrial strength, industry-standard RS485
- EscoBUS communications protocol
- RJ45 ports at rear of each piece of equipment at the rear
- RJ45 to USB interface, plugs into PC USB port
- RJ45 to RS232 interface, plugs into PC DB9/Serial Port
- Each device on the network is assigned with a unique EscoBUS address
- Up to 4 devices of CO₂ Incubator may be daisy-chained to a single PC
- Up to 32 combined devices may be daisy-chained to a single PC
- The maximum distance between equipment is 100 m (328 ft), while the maximum distance between the PC to the last equipment is 1000 m (3280 ft).

System Requirements

Operation system: 32-bit and 64-bit Windows 2000, XP, 7, and Windows 2003 Server

Memory: 512 MB

Hard disk space: 1 GB

Interface : USB/RS232/ RS485

Prerequisite Software: .NET Framework 2.0 (included)

Windows Installer 3.1 (included)

**Internet connection required for email alert feature.*

ORDERING INFORMATION

Software	Item Code	Description
Voyager® and iRMAP™	5250001	Voyager® Software Kit with iRMAP™ subscription for 2 years




iRMAP™

(Intelligent Remote Monitoring Application Protocol)

24/7 Sample Monitoring

Intelligent Remote Monitoring Application Protocol (iRMAP™) is a software system for remote monitoring of Voyager®-compatible Esco devices. A user purchases Voyager®, installs it on a PC and connects it up to the network of Esco equipment. The monitoring data is encrypted and transmitted to an Esco server. User will be able to monitor the freezers/ cabinets 24 hours a day, 7 days a week from anywhere in the world.

The iRMAP, in conjunction with the Predictive Failure Monitoring System can also monitor the health of the Esco unit and proactively dispatch service to correct pending and future issues that could threaten customer samples. iRMAP™ can be accessed anytime, anywhere using a web browser/ iOS/ android mobile application.



PC or Laptop with Voyager® Software



PC or Laptop



Android Phone



iOS Phone
(Beta)




Compatible Esco Equipment

KEY FEATURES

Sample monitoring anytime, anywhere

- 24/7 sample monitoring via internet browser and mobile iOS/ android application
- In the event that your samples are at risk, user will be consistently notified until situation is resolved.
- Every connected Esco unit's health will be monitored and service will be proactively dispatched to correct pending and future issues that could threaten customer samples.



Accessibility


- Users can access Voyager data logs anytime and at any location using either web browser applications, iOS, and android mobile devices.



PC / Laptop Browser




Android Phones



iOS Phones
(Beta)


Alarm

- Alarm is activated if a user-defined parameter limit has exceeded.
- On-screen display of alarm notification and email prompts
- Documentation of all alarms with time and date stamp



Data Log History Display and Graph Functions

- Remote viewing of data log history from multiple devices
- Actual data graph, event and alarm logs history
- Easy export of data as Excel csv file



Voyager® Box

System Overview


Esco Voyager Box is a web-based application developed for remote monitoring, and datalogging of Esco controlled environment laboratory equipment. The device will enable you to monitor your important equipments wirelessly, wireless Local Area Network (LAN) is needed to access the monitored unit. Each monitored parameters can be viewed using web browser of PC, laptop or smartphone. In any event of equipment failure, the application is able to send out remote alarm notification through SMS or e-mail.



Compatible Equipment

- Lexicon II® – Ultra-Low Temperature Freezer (Silver, Gold, and Platinum)
- CelCulture® – CO2 Incubator (CCL)
- CelMate® – CO2 Incubator (CLM)
- Isotherm® – Forced Convection Laboratory Oven
- Isotherm® – Natural Convection Laboratory Incubator (INA)
- Isotherm® – Forced Convection Laboratory Incubator (IFA)
- Isotherm® – Refrigerated Incubator (IFC)





Home Menu

- Touch screen 3,5" display
- Display all monitored parameters
- Access all menu, Wifi connection, add device

Device	
CLS Lab Micro	▲
CCL Lab Electro	▲
IFA Lab Biology	■
IFC Lab	▼
ULTF Lab	▼

Navigation icons: ←, +, ✎, 🗑️

Add Devices


- List monitored device
- Configure measurement
- Add, edit, and delete monitored device

Measurement	
Time	Log
2017-10-13 08:05	Temp: 28.1 / CO2: 5.0 / O2: 20.0 / RH: 60
2017-10-13 08:07	Temp: 27.8 / CO2: 5.0 / O2: 20.1 / RH: 61
2017-10-13 08:09	Temp: 28.1 / CO2: 5.0 / O2: 20.1 / RH: 61
2017-10-13 08:11	Temp: 28.0 / CO2: 5.0 / O2: 20.0 / RH: 60

Navigation icons: ←, ◻, ▶, ☰

Log Measurement

- Measurement history
- Up to 365 days measurement
- Save logging data with USB




Web Interface

- Access all the monitored device in web browser
- Configure the measurement settings as simply as in VoyagerBox device



Installation Requirements



Minimum System Requirements

- Operation system: 32-bit and 64-bit Windows 2000, XP, 7
- Memory: 512 MB
- Hard disk space: 1 GB
- Interface: Wifi Connection

PROtect

A Completely Independent and Redundant Sample Monitoring System

Esco PROtect monitoring system includes radio frequency recorders which can be connected to any brands of freezer, cold storage system and temperature controlled equipment. These recorders automatically send data to a central server to monitor temperature and other parameters and send real-time alert to users in case of severe failure. Esco PROtect complies with quality assurance procedures that is ISO 17025, ISO 15189 norm, GMP and audit trail traceability for conformance to 21 CFR Part 11.

Esco PROtect Wireless Monitoring System can be installed in the following industries/ sectors:

- Hospitals
- Pharmaceutical Industry
- Blood Bank Centers
- Food Safety Laboratories
- Animal Testing Laboratories
- Drugstores
- Research Centers
- Bio-Banking and Repositories
- Private Clinical Laboratories
- Food Industries
- Universities/Academia
- Cosmetic Industry



LEGEND



ESICO Protect
UN/DG



ESICO Protect
ALARM



Visual and
Audible Alert



ESICO Protect
EMD-1



Email Alert

———— Ethernet network

----- Radio Communication
(200-300m)

KEY FEATURES



Real Time Monitoring System

- Wireless Technology that protects your samples 24/7 with continuous data logging function and immediate alarm notifications
- Centralize monitoring for all your equipment (any brands and any models):
 - Temperature for refrigerators, freezers, cold chambers, cryogenic tanks and -80°C freezer
 - Temperature and CO₂ for ovens/incubators
 - Temperature and humidity for climatic chamber, ambient area
- Measurement and recording of physical parameters such as temperature, CO₂, and humidity
- Direct reading on the display
- Radio frequency has longer battery life and better transmission distance, as compared to Wi-Fi
- In using the radio frequency, repeaters can be installed to extend the distance which cannot be easily done by some Wi-Fi system

Simple and Cost Effectiveness

- Fast Installation, commissioning and maintenance
- Can centrally monitor integrity of samples even during a power outage or electronic failure of refrigerator, freezer, incubator, oven, etc
- Accurate and efficient compared to manual data readings
- Reduce labor cost for manual data readings. Manual data readings are less accurate, less reliable and time consuming. Saves you time and money.
- Readily accessible thru the display or monitoring software thru PC
- Multi-site installation
- Interchangeable digital probes with embedded memory to facilitate maintenance procedures.
- Modern and eco-friendly




Traceability and Security – Quality Assurance

- Comply with quality assurance procedures: ISO 17025 and ISO 15189 norms, Good laboratory practices and local regulations recommend the use of tools to guarantee the reliability of your temperature controlled equipment (ultra-low temperature freezers, lab freezers, lab refrigerators, lab ovens, lab incubators, refrigerated incubators, CO₂ incubators, etc).
- Audit trail traceability for conformance to 21 CFR part 11
- Local regulations or recommendations within Quality Management Programs




Alarms

- Real time alarm management
- Warns if problem occurs via audible and visual alarm devices
- Remote alarm notification via SMS/ E-mail




Remote Monitoring and Datalogging

- Remote visualization of data
- Downloaded database is presented in easy-to-understand reports or graphs
- PROtect software enables to visualize recorders located on the installation map and to manage your alarm



System Components



Potential Users of Esco PROtect Wireless Monitoring System

- Hospitals
- Pharmaceutical Industry
- Blood Bank Centers
- Food Safety Laboratories
- Animal Testing Laboratories
- Drugstores
- Hospitals
- Pharmaceutical Industry
- Private Clinical Laboratories
- Food Industries
- Universities/Academia
- Cosmetic Industry

System Requirements for PROtect Software

PROtect Server:

Computer Architecture: 32 bit and 64 bit
 OS: Windows 2000, XP, Vista, Win 7 and Windows 2003-2008 Server
 RAM: 3 GB (min)
 HD: 250 GB capacity or higher
 Video Card: VGA compatible 1280 X 1024 32 bit (minimum)



PROtect Client:

Computer Architecture: 32 bit and 64 bit
 OS: Windows 2000, XP, Vista, Win 7 and Windows 2003-2008 Server
 RAM: 2 GB (minimum)
 HD: 1 GB capacity or higher
 Video Card: VGA compatible 1280 X 1024 32 bit (minimum)

PROtect can be used on the following regions:

- Europe – 868 Mhz
- Africa – 868 Mhz
- Asia – 868 Mhz
- India – 868 Mhz
- ANZ – 902 Mhz
- North America – 902 Mhz
- South America- 868 Mhz/902 Mhz
- Middle East - 868 Mhz

* Currently not applicable for South Korea and Japan region

GENERAL SPECIFICATION		
	Esco PROtect UN	Esco PROtect DG
Measurement Range	From -200 °C to + 400 °C / 4-20mA / 0-1V /dry contact	T : -40 °C to +80 °C / tH : -30 °C to +70 °C from 0 to 100% RH non-condensing
Channel Number	1 or 2	1 or 2
Resolution	0,1°C	0,1°C
Operating Conditions	From -30°C to +70°C	From -30°C to +70°C
Power Supply	Lithium battery	Lithium battery
Battery Lifetime	2 years - Replaceable	2 years - Replaceable
Type of Input	External (Pt 100/Pt 1000/CO ₂ Probe/ 4-20mA/0-20mA/0-1V/dry contact)	External (temperature, relative humidity and dry contact)
Recording Interval	1 second to 90 minutes	1 second to 90 minutes
Internal Memory	20 000 measurements	20 000 measurements
Radio Range	300 m	300 m
Radio Bande	868 Mhz or 902 Mhz ISM	868 Mhz or 902 Mhz ISM
Power	25 mW	25 mW
Digital Display	Yes	Yes
Dimension	123 x 69 x 35 mm	123 x 69 x 35 mm
Protection Index	IP65	IP65
Conformities	EN 12 830, ROHS, CE, FCC, ETS 300-220	EN 12 830, ROHS, CE, FCC, ETS 300-220
Part number	868 Mhz UN-1: 1081050 UN-2: 1081052 902 Mhz UN-1: 1081051 UN-2: 1081053	868 Mhz DG-1: 1081054 DG-2: 1081056 902 Mhz DG-1: 1081055 DG-2: 1081057
Included	User manual, wall-mounting bracket and cable input protection	User manual, wall-mounting bracket and cable input protection

Ordering Information









Modems, Relays and Alarm System

Product Image	Item Codes	Description	Features/Specification
	1081050	Esco PROtect UN-1 868 Mhz band	<ul style="list-style-type: none"> • Single Channel RF Transmitter, -200°C to + 400°C; 4-20mA; 0.1V ; Dry Contact
	1081051	Esco PROtect UN-1 902 Mhz band	
	1081052	Esco PROtect UN-2 868 Mhz band	<ul style="list-style-type: none"> • Dual Channel RF Transmitter, -200°C to + 400°C; 4-20mA; 0.1V ; Dry Contact
	1081053	Esco PROtect UN-2 902 Mhz band	
	1081054	Esco PROtect DG-1 868 Mhz band	<ul style="list-style-type: none"> • Single Channel RF Transmitter, -40°C to +80°C / tH : -30°C to +70°C from 0 to 100% RH non-condensing
	1081055	Esco PROtect DG-1 902 Mhz band	
	1081056	Esco PROtect DG-2 868 Mhz band	<ul style="list-style-type: none"> • Dual Channel RF Transmitter, -40°C to +80°C / tH : -30°C to +70°C from 0 to 100% RH non-condensing
	1081057	Esco PROtect DG-2 902 Mhz band	
	1081058	Esco PROtect MD-1 868 Mhz band	<ul style="list-style-type: none"> • Connects to your PC, it enables communication with Esco PROtect UN/DG • Communicates with an unlimited number of recorders. • Remote alarm driven by integrated dry contact output (vocal telephonic transmitter, buzzer, light). • Software, PC or power supply failure alarm. • Can be used with or without Esco PROtect Relay.
	1081059	Esco PROtect MD-1 902 Mhz band	
	1081060	Esco PROtect EMD-1 868 Mhz band	<ul style="list-style-type: none"> • Enables Ethernet connection to maximize distances between recorders and PC servers
	1081061	Esco PROtect EMD-1 902 Mhz band	
	1081062	Esco PROtect RL-1 868 Mhz band	<ul style="list-style-type: none"> • Enables direct RF communication to maximize distances between recorders and PC servers. • Works with external power supply and supplied with a rechargeable battery. • Includes two power ranges: <ol style="list-style-type: none"> a) 25 mW to communicate with recorders b) 500 mW to communicate with Modem
	1081063	Esco PROtect RL-1 902 Mhz band	
	1081064	Esco PROtect AM-1 868 Mhz band	<ul style="list-style-type: none"> • Provides real time warning if a problem occurs • Integrated sound and visual alarm. • Works with an external power supply and supplied with a backup battery. • Push button for local alarm acknowledgement. • Relay output to trigger remote alarm (vocal telephonic transmitter) • Warns Esco PROtect software when the battery is low.
	1081065	Esco PROtect AM-1 902 Mhz band	

Ordering Information

Options and Accessories			
Product Image	Item Codes	Description	Features/Specification
	1081066	PROtect Test-1 868 MHz	<ul style="list-style-type: none"> Use as a recorder to test connection between modem and recorder with actual percentage signal
	1081067	PROtect Test-1 902 MHz	
	1081068	PROtect Modem Test-1 868 MHz	<ul style="list-style-type: none"> Use as a modem to test connection between modem and recorder with actual percentage signal
	1081069	PROtect Modem Test-1 902 MHz	
	1081070	PROtect Software and Dongle	<ul style="list-style-type: none"> Remote visualization of data. Downloaded database is presented in easy to understand reports or graphs. PROtect software enables to visualize recorders located on the installation map and to manage your alarm. Compliant with FDA (CFR 21 part 11) <p><i>Note: PROtect software and dongle needs to be installed even if one (1) PROtect recorder (UN or DG) is connected to the equipment.</i></p>
	1081071	Audible and Visual Alarm Device	<ul style="list-style-type: none"> Warns if an alarm occurs. Works with an external power supply. Connects directly to PROtect Modem or to PROtect Alarm. Audible and visual Alarm.
	1081072	PT100 Probe Class B	<ul style="list-style-type: none"> Measurement range: -200°C to + 260°C Type of cable: 3 wires (3m)
	1081288	PT100 Probe Class A	
	1081289	PT100 Probe Class 1/3 DIN	
	1081073	Digital Temperature Probe	<ul style="list-style-type: none"> Type of probe: Dallas Measurement range: -40°C to + 80°C (-45°C max) accuracy: ± 0.3°C Type of cable: flat 3m cable Embedded memory for gauging parameters storage
	1081075	Battery for Recorders	<ul style="list-style-type: none"> Lithium Battery for 10 pcs
	1081170	Digital Temperature / Humidity Probe	<ul style="list-style-type: none"> Measurement range : -40 +80°C /0 100%HR Type of probe for T°: Pt 100 Class A Class A Accuracy: ±2%HR from 0 to 90%HR ±3%HR from 90 to 100%HR Power: delivered by the SPYRF recorder
	1081169	Digital Temperature/ Humidity Probe with 3m Extension Lead	
	1081293	Digital Temperature/ Humidity Probe with 12m Extension Lead	
	1081166	Digital (Low) Temperature	<ul style="list-style-type: none"> Type of probe: Low Temp Stainless Steel probe Type of cable: 3m Teflon Cable Measurement range from 0°C to -200°C accuracy: ± 0.2°C Immersion of cable: continuous immersion in water and glycolic solutions and punctual immersion in alcoholic solutions"
	1081080	Universal Adaptor Cable for CO ₂ and DP sensors	<ul style="list-style-type: none"> (5m) with connector IP67 for PROtect UN and DG.

Ordering Information


Options and Accessories			
Product Image	Item Codes	Description	Features/Specification
	1081165	Extension Lead for Probe	<ul style="list-style-type: none"> • 3m Extension Lead for Digital Temperature Probe
	1081262	SMS Modem	<ul style="list-style-type: none"> • Warns if any alarm occurs by sending SMS to receiver's number • Needs a SIM card • Supports multiple receiver's number • USB Connection
	1081076	PROtect CO ₂ Probe	<ul style="list-style-type: none"> • Measurement range : 0 – 20% • Output : 4-20 mA • Mounting : Wall mounting bracket supplied • Supply: 100-240 VAC
	5072029	PROtect Differential Pressure Probe	<ul style="list-style-type: none"> • Measurement Range: -100 Pa to +100 Pa • Output: 4 - 20 mA • Supply: 100-240 VAC (include Power Supply 1081263)
	1081078	Power Supply for Modem	<ul style="list-style-type: none"> • Power Supply for MD-1, EMD-1, AM-1
	1081079	Power Supply for Relay	<ul style="list-style-type: none"> • Power supply for RL-1
	1081167	PROtect Door Contact (With Sliding Doors)	<ul style="list-style-type: none"> • Use to sense door contact that door is open or close
	1081168	PROtect Door Contact (Non-Sliding Doors)	<ul style="list-style-type: none"> • Use to sense door contact that door is open or close
	5072009	PROtect GSM Auto Dialer set	<ul style="list-style-type: none"> • Warns if any alarm occurs by call receiver number (auto dialer) • Needs a SIM card • Supports until three (3) receiver numbers • Supply: 100-240 VAC (include Power Supply 1081291)

SampleTracker

System Overview

SampleTracker is a comprehensive Bio sample Inventory Management solution. The Software tracks the sample management activities throughout its life cycle. The software scope covers all the features and functionalities required for the Sample Inventory Management and Tracking for all sizes of laboratories, managing different types of samples and having various types of storages. Flexibility in configuration of laboratories and sample attributes is one of the unique characteristics of SampleTracker.

SampleTracker automatically generates a unique identification number for each sample, which makes it easy for users to trace the sample as and when required



KEY FEATURES

Storage Template

Select a storage type to define storage template. Specify storage layout and provide details for the template.

Storage Type: Storage Template Name: Storage Template Type: Activate

Enforce Barcode

Storage Capacity

Columns: Rows: = 64 Number of storage

Specify number of rows and columns to define storage capacity:

Std Box 64							
1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40

Storage Details

Other details

Manufactured by:

Box Properties

Color: Purchased on:

Box Details

Box name: Capacity:

Laboratory Storage Management

- Detail storage definition and configuration
- User able to manage their sample such as save, update and save as another storage
- User will allow to design their own storage template such as capacity and details storage information

Configuration Template

Delete Template

Esco Freezer configuration

- Esco freezer 1
 - Compartment 1
 - Rack 1
 - Rack 2
 - Compartment 2
 - Compartment 3
 - Compartment 4

Template Details

Storage Overview View Details

Rack 1 (Rack, 4*5)

Storage Capacity -20
Number of cells equal to the storage capacity.

1	2	3	4	5
6	7	8	9	10
11	12	13	14	15
16	17	18	19	20

Lab Management

- Detail hierarchies as a configuration template
- navigate through all the configuration templates listed

New Query

Create a query using the given options. Click Run Query to see the results in result table.

Name: Query Type: Query Focus:

And/Or	Field	Operator	Field Value
<input type="button" value="Click to add clause to query"/>			

The results can be seen in Storage Management view to take further actions.

Query Results

Total Results - 0 | | | |

Query Builder

- User will allow to create their own query to search the specific data from database
- User can save the most popular queries and can use these queries in future instead of creating new query every time



User allow to mark their sample that will help them for identifying the individual label in the million samples



SampleTracker / Alto Platinum Controller functions



Login to SampleTracker

At the Alto Platinum controller home screen go to application and select Inventory System

Type in username and password at SampleTracker login screen.



Performing SampleTracker Functions

- **CHECK IN:** It allows user to Check In the samples again into the system after processing. Only Checked Out samples can be Checked In.
- **CHECK OUT:** It allows user to take out the samples for processing. Only Checked In samples can be Checked out.
- **DISPOSE:** Dispose option deletes the selected samples. This option can be used to remove unwanted or degraded samples so that these positions are available for new samples. Any one of Checked In, Checked Out or Moved Out State allows disposal of the samples.

Minimum Requirements

- .Net Framework 3.5 to 4.0
- SQL Server 2008 R2
- Windows 7 SP1, SP2
- Hard Drive 80 – 100 GB
- RAM 1GB

After Sales Services



Parts Availability

Whenever service is needed and parts are required, minimizing downtime is a critical objective. Statistical usage analysis helps Esco to predict parts life, permitting Esco to manage logistics and stage proper inventories around the world. The combination of predictive maintenance, historical data and geospecific proximity assures our customers that parts and labor are available whenever service is scheduled through the local sales organization.

Registration, Documentation and Instruction

Quality control at Esco extends from research and development through engineering, manufacturing, shipment, delivery and customer feedback. Esco maintains an aggressive program to encourage warranty card registration by mail, email or online submission so that we know where Esco products are located and how they are being used. Rest assured that all information disclosed from warranty registrations will be kept confidential. All Esco products include unique serial numbers for identification. Documentation for all performance tests is archived and maintained for customer reference.

Online Technical Information

Site preparation instructions are useful before product arrival and installation. Installation and start-up manuals, operation manuals and quick reference guides are available anytime from the Esco resources online. An interactive online LiveSupport™ concierge center accessible through the Esco website offers extended hours of operation. LiveSupport™ permits users to dialogue directly with Esco personnel.

NSF International Accreditations

The National Sanitation Foundation (NSF) International is an independent, non-profit organization that provides standards development, product certification, auditing, education and risk management for public health and the environment.

In line with Esco's commitment in providing world class services worldwide, Esco has a large contingent of NSF accredited certifiers which makes Esco not only an **Excellent Standards COmpany** but also an **Excellent Service COmpany**, which exemplifies Esco's collective quest of being an **Eternally Successful COmpany**.

The NSF mark is your assurance that the product complies with all the standard requirements, tested by one of the most respected independent certification organizations in existence today. NSF conducts periodic unannounced inspections and product testing to verify that the product continues to comply with the standard. It is valued by consumers, manufacturers, retailers and regulatory agencies worldwide.


References and Links

For more information, you can visit Esco at www.escoglobal.com

Whatever your needs are,
**we serve it to you
on a silver platter.**



ESCO GLOBAL NETWORK



- 📍 Global Offices
- 📍 Joint Ventures
- 📍 Licensee
- Distributors
- Factories
- R&D Centers
- Regional Distribution Centers



- ART Equipment
- Biological Safety Cabinets
- CO₂ Incubators
- Compounding Pharmacy Equipment
- Containment / Pharma Products
- Ductless Fume Hoods
- Freeze Dryers
- Lab Animal Research Products
- Laboratory Centrifuges
- Laboratory Fume Hoods
- Laboratory Ovens and Incubators
- Laboratory Shakers
- Laminar Flow Clean Benches
- PCR Cabinets
- PCR Thermal Cyclers
- Powder Weighing Balance Enclosures
- Ultra-low Temperature Freezers

The Esco Group of Companies is a global life sciences tools provider with sales in over 100 countries. The group is active in lab equipment, pharma equipment and medical devices. Manufacturing facilities are located in Asia and Europe. R&D is conducted worldwide spanning the US, Europe and Asia. Sales, service and marketing subsidiaries are located in 12 major markets including the US, UK, Singapore, Japan, China and India. Regional distribution centers are located in the US, UK, and Singapore.

Life Science • Chemical Research • Assisted Reproductive Technology (ART) • Pharmaceutical Equipment • General Equipment



WORLD CLASS. WORLDWIDE.

Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777
 Tel +65 6542 0833 • Fax +65 6542 6920 • mail@escoglobal.com
 www.escoglobal.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA
 Toll-Free USA and Canada 1-877-479-3726 • Tel 215-441-9661 • Fax 484-698-7757
 eti.sales@escoglobal.com • www.escolifesciences.us

Esco Global Offices: Bangladesh | China | Denmark | Hong Kong | India | Indonesia | Italy | Japan | Lithuania
 Malaysia | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam

